

myLab Modules (New Web Version) – **MUHC** RAD User Overview

Module	Menu item (if applicable)	Used for...
Inventory	Quick Isotope Entry	<ul style="list-style-type: none"> ▪ Adding new radioactive materials to your inventory
	Disposal/Add Usage	<ul style="list-style-type: none"> ▪ Recording the usage of radioactive materials in your inventory ▪ Transferring radioactive materials from your inventory to another PI or to record aliquots ▪ Editing or deleting a usage/disposal previously recorded
	Package Receipt Survey NEW	<ul style="list-style-type: none"> ▪ Recording results of performed wipe tests upon receipt of items
	View In-Lab Waste Containers	<ul style="list-style-type: none"> ▪ Preparing a container for waste pick-up – closing a container + recording results of a performed wipe test (NEW) + selecting a temporary waste area where the container will be dropped off
PI Overview NEW		<ul style="list-style-type: none"> ▪ Viewing information about permits, training, lab/room maps (<i>if available</i>), and inspection results

User Guides are available in the [IT Knowledge Base](#).